



CONTACT US:

Office :

WenZhou YouMai International Trading CO.,Ltd.

Address:Room A 1702-1,Xinyi Building,Shifu Road,Wenzhou, Zhejiang, China

Tel:+86-577-88267991

Web:www.unimaster.com.cn

Factory :

Address:No.1568,Binhai 1st Road,Wenzhou Economic and Technology Development Zone

Tel:+86-577-86805999

Cannula



- Sharp RS type



- Cannula U type



- Cannula L type



- Cannula W type



- Taper Point (Blunt)

UNIMASTER[®] THREADS COLLECTION

Unimaster offer various types of lifting threads for a number of different lifting procedures on the body and face.

The threads range includes:
The UP, Embossul Sawtooth, The Fishbone, U Cog, The Nose, The Lift, The Multiple, The Screw, The Mone, The Noon, Screw Tornado

Material-wise threads are available in PDO, PCL


Thread Type	PDO	PCL
Biodegradation period	6-9 months	6-9 months
Duration of effect	Up to 12 months	18-24 months

Mesothreads

Mesothreads are available in 2 different materials: PDO, PCL.

All the materials have been widely used in medicine for decades and are safe for aesthetic medicine procedures.

- [1.The Up](#)
- [2.Embossul Sawtooth](#)
- [3.The Fishbone](#)
- [4.U Cog](#)
- [5.The Nose](#)
- [6.The Lift](#)
- [7.The Multiple](#)
- [8.The Screw](#)
- [9.The Mono](#)
- [10.The Noon](#)
- [11.Screw Tornado](#)



Sharp RS type · Cannula U type · L type · **W type** · Taper Point (Blunt)

The Up

Lift the sagging face, buttocks, chest. W-shaped needle, 18-19G, USP 2 / 0, rotary single barbed serrated.

Type	Material	Needle size	Thread length	USP	Cannula
The Up	PDO/PCL	18 G 100mm	150 mm/160mm	2	W Type
	PDO/PCL	19 G 100mm	150 mm/160mm	0	W Type



Embossul Sawtooth

Longer-lasting lifting effect. W-shaped needle, 18-19G, USP 2, single side without spiral embossing serrated thread.


Type	Material	Needle size	Thread length	USP	Cannula
Embossul Sawtooth	PDO/PCL	18 G 100mm	160 mm	2	W type
	PDO/PCL	19 G 100mm	160 mm	2	W type



The Fishbone

A longer-lasting lifting effect. W-shaped needle, 18-19G, USP 2, the two ends of the herringbone line are in opposite directions and there is a 90° turn in the middle.


Type	Material	Needle size	Thread length	USP	Cannula
The Fishbone	PDO	18 G 100mm	160 mm	2	W type
	PDO	19 G 100mm	160 mm	2	W type



U Cog

It is used for large-scale lifting of breasts, buttocks, forehead, eyebrows and other parts. Solid needle, 19-20G, USP 2 / 1, Swivel Serrated Wire.


Type	Material	Needle size	Thread length	USP	Cannula
U Cog	PDO	19 G 140mm	410 mm	2	Taper point
	PDO	20 G 100mm	410 mm	1	Taper point



The Nose

Lift the nose (Nose Bridge/Nose Column). L-shaped needle, 19G, USP 2, single thorn serration without rotation.


Type	Material	Needle size	Thread length	USP	Cannula
The Nose	PDO/PCL	19 G 38mm	50 mm	2	L type
	PDO/PCL	19 G 60mm	80mm	2	L type



The Lift

Improve facial sagging, lift mid-face, lift eyebrows, lift eye corners and other facial local enhancements. L-shaped needle, 21G, USP 2-0, single-barb serration.

Type	Material	Needle size	Thread length	USP	Cannula
The Lift	PDO	21G 60mm (single Line)	120 mm	2-0	L type
	PDO	21G 60mm (Double Line)	200 mm	2-0	L type



The Multiple

Mono threads are singular, smooth sutures that combat ptosis and stimulate neocollagenesis processes by creating an intradermal framework.

Type	Material	Needle size	Thread length	USP	Cannula
The Multiple	PDO/PCL	23 G 38mm	50 mm	5-0	L type
	PDO/PCL	23 G 60mm	80 mm	5-0	L type



The Screw

For face and body, jawline, nasolabial folds, mouth lines, fine lines around the eyes/butterfly sleeves, abdomen, buttocks. RS-shaped needle, 26-30G (no 28G), USP 5-0/6-0, helical barbless wire.

Type	Material	Needle size	Thread length	USP	Cannula
The Screw	PDO	27 G 50mm	70 mm	5-0	RS type
	PDO	29 G 25mm	30 mm	6-0	RS type
	PDO	29 G 38mm	50 mm	6-0	RS type
	PDO	30 G 25mm	30 mm	6-0	RS type
	PDO	26 G 60mm	80 mm	5-0	RS type



The Mono

Tighten facial tissues and relax muscles. Forehead collapse, apple muscle collapse, temple collapse, facial depression and other partial facial collapse problems/wrinkle removal for forehead lines, Chevron pattern, neck lines and other body parts. RS-shaped needle, 26-30G (no 28G), USP 5-0/6-0, helical barbless wire.


Type	Material	Needle size	Thread length	USP	Cannula
The Mono	PDO/PCL	27G 50mm	70 mm	5-0	RS type
	PDO/PCL	29G 25mm	30 mm	6-0	RS type
	PDO/PCL	29G 38mm	50 mm	6-0	RS type
	PDO/PCL	30G 25mm	30 mm	6-0	RS type
	PDO/PCL	26G 60mm	80 mm	5-0	RS type



The Noon

Wrinkles around the eyes, loose bags under the eyes, tear grooves and other aging problems in the eye area. U-shaped needle, 30G, USP 6-0, spiral without thorn thread.

Type	Material	Needle size	Thread length	USP	Cannula
The Noon	PDO/PCL	30G 25mm	30 mm	6-0	U type



Screw Tornado

Corrects wrinkles in areas with variable facial expressions. Spiral double wire, RS-shaped needle, 26-27G, USP 7-0, two spiral wires.

Type	Material	Needle size	Thread length	USP	Cannula
Screw	PDO	26 G 60mm	90 mm	7-0	RS type
Tornado	PDO	27 G 50mm	70 mm	7-0	RS type



Paving The Safe Way For Aesthetic