

**unimaster<sup>®</sup>**  
PAVING THE SAFE WAY FOR AESTHETIC

**CONTACT US:**

Address: Room A 1702-1, Xinyi Building, Shifu Road, Wenzhou, Zhejiang, China

Email: [Info@Unimaster.com.cn](mailto:Info@Unimaster.com.cn)

[www.unimaster.com.cn](http://www.unimaster.com.cn)

# Medical Cosmetic Series

**unimaster<sup>®</sup>**





## About Us

Driven by a sense of “ Deliver Best Care to the World”, for nearly half-century, we have been enthusiastic and unwavering on the journey to safer injection and better care for medical.

From 1980’ s, with the first piece of medical disposable injection needle was created in Wenzhou, we started to involve and became the industry pioneer in needle manufacturing. And we delivered the innovation and contributed to move the medical injection field of China forward to be safer and more efficient.

We have leveraged deep expertise and advanced technologies to continuously develop in our needles and since 2015, we have been the first company in China, who entered into the research on aesthetic injection and succeeded in launching the Nano Needle technology and treatment conception for medical beauty applications. Now our products are diligently designed and diversified to enable different people who are pursuing of the beauty to benefit from the reliable and effective drug delivery.

Unimaster, an ethic expression for faith, endows us the DNA in advancing to empower more and more people to have an expected living and unlimited sense of being. We see more than your beauty, We see you !



# Why choose Unimaster



With strict verification system, we select premium Japan Stainless Steel as the fundamental to ensure the high quality of our needles with maximum sharpness and less corrosion.



Our unique and advanced Nitrogen Hardening Technology in handling on the bevel to realize super performance in 20, 30, 50 multi times injection without damage to needle tip. Maximum 200 times puncture.



30 Years know-how from strip welding, extruding, OD Reduction, Cutting, Polishing, Grinding and Ultrasonic Washing on cannula to coated and assembled needle, our technology allows for smoother insertion into the skin.



Other than TW and RW conventional needles, we are able to make Ultra-Thin Wall for nano needles to create maximum lumen and improve the flow rate of liquid by more than 20%.



# Meso Needle



## 1.2mm Ultra Short Needle

- Unimaster needle is made of hard Japanese steel that resists blunting.
- The micro-polishing technology offers a superior sharpness and ultra-thin wall structure.
- It maximizes patient comfort with an almost painless injection experience while minimizing the effects of bleeding, bruising, swelling.



- Adjustable series
- Adjust the length range to 0-7mm by rotating the screw knob.

### Specification:

Metric Size(mm)	0.3	0.26	0.23	0.2	0.18
Inch Size	30G	31G	32G	33G	34G
Hub Color	Yellow	White	Green	Grey	Orange
Nominal Length	1.2mm-25mm				
Wall Thickness	Ultra Thin Wall				
Packing	100pcs/box, 50boxes/carton, 5000pcs/carton				

## Multi style packaging (recommended) :



- Plastic package
- High cost performance  
High market share

- Premium package
- No paper scraps  
No risk of lesions of vessels

- Tasse package
- Innovative, Unique  
More easier stripping




# Micro Needle



Anatomical Region		Cannula recommendation	
Periorbital		25, 27, 30G	25 ~ 40mm
Forehead		25, 27G	40 ~ 50mm
Brow		25, 27G	25 ~ 50mm
Temporal, Middle third of the face, Lower third of the face		25, 27G	40 ~ 50mm
Lips		25 ~ 30G	25 ~ 40mm
Chin, Neck, Decollete, Hands		25 ~ 27G	40 ~ 50mm

17G	70 / 90mm
18G	50 / 60 / 70 / 90mm
19G	50 / 60 / 70 / 90mm
20G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
21G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
22G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
23G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
24G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
25G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
26G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
27G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
29G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
30G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm

## Detail display






# Multi needle — 4 Pins

## Water Mesotherapy Working Theory

Water mesotherapy is to inject hyaluronic acid into skin. Hyaluronic acid contains moisture,200-300 times as rich as one's body. So that it makes the skin smooth, soft, bright and hydrated after the injection.




## Specification


Metric Size(mm)	0.23	0.18
Inch Size	32G	34G
Hub Color	Transparent color	
Nominal Length	1.0mm~1.5mm	
Wall Thickness	Ultra-Thin wall	
Packaging	20pcs/box, 60boxes/carton,1200pcs/carton	

Note: Models are classified by product packaging


## Detail display



## Effect display:



## Principle of action:





# Multi Needle

## INSTRUCTIONS:

1. According different treatment area and ficial cases to choose different multi pins.
2. Beauty products will not leak due to negative pressure. The best choice for skin care treatment
3. Sharp needle tips gives smaller wounds and speeds up recovery. Improve wound healing and anti-aging ability.
4. One skin booster for multiple use, reduces number of times needed for puncture and pain

## Dermal Fillers

## Treatment Areas



## SERIES:

### Dermal Series

Machine 16/9/5Pins

Adjustable 16/9/5Pins

Visual dial 16Pins

### Crystal Series


Machine 16/9/5Pins

Adjustable 16Pins

Hand 16/9/5Pins

• Transparent protective sleeve

• Sharp pinhole



• visual dial,  
The operator is more clear about the exposed length,  
More accurate delivery of the required dose.






# PREFILLED SALINE FLUSH SYRINGE (0.9% NaCl)

## Product Introduction:

- Free of latex, DEHP, PVC, BPA and Preservatives
- Medical Grade Polypropylene
- USP-qualified Saline Solution
- Brominated Butyl Rubber for Stable Performance
- Termimnally Sterilized with Steam-heat



## Application

Prefilled Saline Flulsh Syringe (0.9% NaCl) is intendedand solely intended for flushing and locking indwelling vascular access devices after medications are given.It is compatible with inwelling needles, PICCs, CVCs and infusion ports.It is not intended for infusion.

## Advantages

Helps Reduce Preparation Time, Increase Efficiency

Reduces Bacterial Reproduction and Risk of Infection

Reduces Particulate Contamination

Prevents Blood Clots

Reduces medicine interactions between various medications

Moderate Flushing Pressure Lowers Damage to Catheter Wall

## Specification

Reference Code	Item	Size
HP-FO01	Silicone Oil, PFS	N/A
HP-FO04	Barrel, PFS	10ml
HP-FO05	Tip Cap, PFS	N/A
HP-FO06	Bromobutyl Rubber Plunger, PFSa	N/A
HP-FO07	Polypropylene Plunger Rod, PFS	N/A
HP-C3B002	0.9% NaCl, PFS	10ml
HP-C3B004	0.9% NaCl, PFS	5ml



# Needle Holder



■ Product Introduction:

- Metal material
- Easier to unscrew needle
- Multiple use
- Match all size needle

# Cosmetic Syringe

■ Advantages:

- PC raw material
- Higher Quality than normal syringe
- Match different liquid concentration
- Low dead space design



Wider finger flange for easier puncture  
Crystal plunge clearly to see the fluid



We Always Make Great Effort  
For  
Human Health



**unimaster**<sup>®</sup>

WENZHOU  
YOU MAI  
INTERNATIONAL TRADING CO.,LTD.